

GCE HEALTHCARE

Zen-O™

ENJOY FREEDOM AND INDEPENDENCE
– WITH MOBILE OXYGEN ON DEMAND.



GCE[®]
Healthcare

EASY MAINTENANCE

Zen-O™ is incredibly easy to maintain. Our design makes it simple for your home oxygen provider or retailer to carry out basic maintenance quickly and easily.

BUILT TO LAST

Zen-O™ meets the requirements of the European Medical Device Directive and the Food and Drug Administration (FDA) in the US.

Zen-O™ comes with a product warranty of three years or 15,000 hours.



OFFERING FREEDOM AND INDEPENDENCE

We understand that when receiving long-term oxygen therapy, it helps to keep an active lifestyle. Zen-O™ is a portable oxygen concentrator that, unlike conventional home oxygen equipment, allows you the freedom to move and enjoy activities within and outside of your home.



RESPONSIVE OXYGEN DELIVERY

In pulse mode, Zen-O™ has best in class breath detection technology, meaning it is able to detect shallow breaths. It also uses our Rate Responsive Therapy (RRT) technology to deliver a constant bolus of oxygen in response to your breath rate. Zen-O™ provides optimum therapy by delivering the oxygen bolus within the first half a second of your breath.



KEEPING YOU ON THE MOVE

Zen-O™ comes with everything you need to remain active, including a discreet carry bag, a pull cart, an accessory bag, AC & DC power supply packs and a re-chargeable lithium-ion battery. Zen-O™ can hold up to two batteries that can last up to an impressive eight hours in pulse mode. Zen-O™ is even suitable for air travel, having met the requirements of the US Federal Aviation Administration (FAA).

DESIGNED FOR EVERYDAY USE

With its intuitive buttons and clear LCD display, you'll find the portable Zen-O™ is simple to operate wherever you are and whatever activity you're doing. Zen-O™ weighs just 4.66kg (10.25lbs), delivers up to 2 litres of oxygen a minute and it offers versatility with two flow options – pulse and continuous.

Zen-O™ is complete with audio and visual alarms that are easy to recognise, providing you with assurance on the performance of your device.





GET IN TOUCH

Would you like to know more? Please visit our website, www.gcehealthcare.com/zen-o, where there is a wealth of information about Zen-O™, as well as our smaller and quieter device Zen-O lite™.