

GCE HEALTHCARE

# Zen-O™

ENJOY FREEDOM AND INDEPENDENCE  
– WITH MOBILE OXYGEN ON DEMAND.



**GCE**<sup>®</sup>  
Healthcare

## EASY MAINTENANCE

Zen-O™ is incredibly easy to maintain. Our design makes it simple for your home oxygen provider or retailer to carry out basic maintenance quickly and easily.

## BUILT TO LAST

Zen-O™ meets the requirements of the European Medical Device Directive and the Food and Drug Administration (FDA) in the US.

Zen-O™ comes with a 3 year or 15,000 hours device warranty.\* Warranty of the sieve bed assembly supplied in the device is limited to 2 years or 15,000 hours.



## OFFERING FREEDOM AND INDEPENDENCE

We understand that when receiving long-term oxygen therapy, it helps to keep an active lifestyle. Zen-O™ is a portable oxygen concentrator that, unlike conventional home oxygen equipment, allows you the freedom to move and enjoy activities within and outside of your home.



## RESPONSIVE OXYGEN DELIVERY

In pulse mode, Zen-O™ has best in class breath detection technology, meaning it is able to detect shallow breaths. It also uses our Rate Responsive Therapy (RRT) technology to deliver a constant bolus of oxygen in response to your breath rate. Zen-O™ provides optimum therapy by delivering the oxygen bolus within the first half a second of your breath.



## KEEPING YOU ON THE MOVE

Zen-O™ comes with everything you need to remain active, including a discreet carry bag, a pull cart, an accessory bag, AC & DC power supply packs and a re-chargeable lithium-ion battery. Zen-O™ can hold up to two batteries that can last up to an impressive nine hours in pulse mode. Zen-O™ is even suitable for air travel, having met the requirements of the US Federal Aviation Administration (FAA).

## DESIGNED FOR EVERYDAY USE

With its intuitive buttons and clear LCD display, you'll find the portable Zen-O™ is simple to operate wherever you are and whatever activity you're doing. Zen-O™ weighs just 4.66kg (10.25lbs), delivers up to 2 litres of oxygen a minute and it offers versatility with two flow options – pulse and continuous.

Zen-O™ is complete with audio and visual alarms that are easy to recognise, providing you with assurance on the performance of your device.





## GET IN TOUCH

Would you like to know more? Please visit our website, [www.gcehealthcare.com/zen-o](http://www.gcehealthcare.com/zen-o), where there is a wealth of information about Zen-O™, as well as our smaller and quieter device Zen-O lite™.

Doc No.: DL-00521 Rev.3 © GCE Healthcare.  
Part of GCE Group, a world leader in gas control equipment.

\* Refer to GCE limited warranty policy for more details